

Kinematic Improvements 1220 N. Coit Road Suite 101 Plano, Tx 75075

Bounders
Free Range
AFO

Patient NameSexHt WtDxCustom name within padding – Yes or No Left Right Bilateral	
Cast Date Need by Date Facility Name Orthotist Address	
CityStateZIPPO Shipping AddressStateZIP	P#
Adjustable Bounder - Posterior Dynamic Bounder Strap –	
 Free Range Bounder comes with a plantar flexion assist strap on the posterior aspect of the AFO. It has 3 slots on the back to change the plantar flexion tension with no tools needed. We will select the correct posterior strap based off the height of the cast and weight of the patient Bounder night time stretching straps YES NO Cast Corrections in Degrees Neutral Dorsiflexion Plantar Flexion Leave as is Cast Modifications to address Pronation Supination Neutral Leave as is 	
	StrappingColor of Straps
Plastic and Pattern Selection Transfer Paper Color Polypropylene Modified-Polyethylene Copolymer Other Thickness of Plastic 1/16" 3/32" 1/8" 5/32" 3/16" 1/4" Anterior Tibial Shell (2/3 rd the length of the tibia) Calf paddingYes No	Proximal Velcro Strap 1" 1 ½" 2" Distal Ankle Strap 1" 1 ½" Attach Strap Lateral Side Medial side Type of closure Chafe Velcro Hook Instep strap – Attach Medial Lateral None Forefoot Strap – Attach Medial Lateral None Figure of 8 Strap
** Perf (standard) Aliplast	Shipping UPS Ground
Padding / Inner Boot Color Selection UPS Next Day Saver UPS Next Day by 10am UPS 2 nd Day Air	

Fax orders to 972-767-3216 or Email them to <u>Kinematicimprovements@gmail.com</u>

Send Negative Cast to Address Above

Call **888-338-8309** for technical questions & order Help